




Oxivir (AHP) Glove Substitute Selection Reference

Product	Maximum Duration
Medline Sensicare ICE 	20 minutes
Medline Sensicare 	20 minutes
Aurelia Sonic 	5 minutes
Medline Sensicare Silk 	10 minutes
VersaShield Extended Cuff** 	35 minutes
Curad Extended Cuff** 	10 minutes


*Gloves are not required for Oxivir (accelerated hydrogen peroxide) use less than 5 minutes. Hand rinsing is recommended.

**Extended cuff gloves should be reserved for tasks specifically identified in the [Non Sterile Exam Glove Selection Guide](#).

Refer to the [Product Change Notice](#) or [Non Sterile Exam Glove Selection Guide](#) for more detailed information and wear times with specific chemicals or hazardous drugs.

Availability subject to change. Current as of October 6, 2020

Oxivir (AHP) Glove Substitute Selection Reference

Product	Maximum Duration
Medline Sensicare ICE 	20 minutes
Medline Sensicare 	20 minutes
Aurelia Sonic 	5 minutes
Medline Sensicare Silk 	10 minutes
VersaShield Extended Cuff** 	35 minutes
Curad Extended Cuff** 	10 minutes

*Gloves are not required for Oxivir (accelerated hydrogen peroxide) use less than 5 minutes. Hand rinsing is recommended.

**Extended cuff gloves should be reserved for tasks specifically identified in the [Non Sterile Exam Glove Selection Guide](#).

Refer to the [Product Change Notice](#) or [Non Sterile Exam Glove Selection Guide](#) for more detailed information and wear times with specific chemicals or hazardous drugs.

Availability subject to change. Current as of October 6, 2020